

Stearic acid-PEG-Mal (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Stearic acid-PEG-Mal (MW 5000) is used for surface modification of various substrates through the reaction between hydroxyl and ethoxysilane groups. This polymer is also suitable for click chemistry reactions with alkynyl-containing molecules. It is reagent grade and intended solely for research purposes.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481