

DOPE-PEG-Amine (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOPE-PEG-Amine (MW 3400) is a phospholipid known for its excellent biocompatibility and distinct amphipathic properties. These characteristics make it an ideal and primary excipient or dosage form in pharmaceutical formulations, enabling enhanced therapeutic outcomes. Additionally, DOPE-PEG-Amine (MW 3400) can be utilized in drug delivery research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481