

DOPE-PEG-Amine (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOPE-PEG-Amine (MW 5000) is a phospholipid known for its excellent biocompatibility and significant amphiphilic properties, making it a key and suitable excipient in drug formulations to achieve enhanced therapeutic outcomes. It is utilized in research focusing on drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481