

DOPE-PEG-Mal (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOPE-PEG-Mal (MW 5000) is a phospholipid known for its excellent biocompatibility and notable amphiphilic properties, making it an ideal and primary vehicle or excipient in drug formulations for enhanced therapeutic outcomes. It is applicable in research related to drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481