

DSPE-PEG-NHS, MW 1000

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG-NHS, with a molecular weight of 1,000, is an amine-reactive phospholipid. The reaction between the NHS ester and an amine results in a stable amide bond. This polymer is amphiphilic, enabling it to form micelles or lipid bilayers in aqueous solutions.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481