

Polyethylene glycol distearate (Mn 930)

Chemical Properties

CAS No. : 9005-08-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Poly(ethylene glycol) distearate (Mn 930) is a derivative of polyethylene glycol, utilized as a crosslinker and in click chemistry reactions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481