

Antibody buffer (20mM HEPES, 150 mM NaCl, 240mM Sucrose, pH 7.2 buffer)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at 4°C  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Antibody buffer (20mM HEPES, 150 mM NaCl, 240mM Sucrose, pH 7.2 buffer) is a commonly used antibody storage buffer for biochemical experiments.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481