

Potassium phosphate buffer (0.1 M, pH 6.9)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: 4°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Potassium phosphate buffer (0.1 M, pH 6.9) is a phosphate buffer solution widely employed in life science research and experimentation.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481