

Poly(ethylene glycol) diacrylamide (MW 2000)

Chemical Properties

CAS No. : 160556-48-9

Formula: (C₂H₄O)_nC₁₀H₁₆N₂O₃)

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Poly(ethylene glycol) diacrylamide (MW 2000) is a reagent used in biochemical reactions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481