

## Methoxypolyethylene glycol succinimidyl succinate

### Chemical Properties

CAS No. : 78274-32-5

Formula: (C<sub>2</sub>H<sub>4</sub>O)<sub>n</sub>C<sub>9</sub>H<sub>11</sub>N<sub>06</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Methoxypolyethylene glycol succinimidyl succinate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481