

## 3'-NH2-ddCTP sodium

### Chemical Properties

CAS No. :

Formula: C<sub>9</sub>H<sub>17</sub>N<sub>4</sub>O<sub>12</sub>P<sub>3</sub>.xNa

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	3'-NH <sub>2</sub> -ddCTP sodium is a labeled and modified deoxynucleoside triphosphate (dNTP) that acts as chain-elongating inhibitor of DNA polymerase for DNA sequencing.
Targets(IC50)	DNA/RNA Synthesis

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481