

dTTP α S sodium

Chemical Properties

CAS No. : 178603-75-3

Formula: C₁₀H₁₇N₂O₁₃P₃S.xNa

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	dTTP α S sodium is a labeled and modified deoxynucleoside triphosphate (dNTP) that acts as chain-elongating inhibitor of DNA polymerase for DNA sequencing.
Targets(IC50)	DNA/RNA Synthesis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481