


(S)-Quinuclidin-3-amine (dihydrochloride)

Chemical Properties

CAS No. : 119904-90-4
 Formula: C₇H₁₆Cl₂N₂
 Molecular Weight: 199.12
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-Quinuclidin-3-amine (dihydrochloride) (Standard) is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0221 mL	25.1105 mL	50.221 mL
5 mM	1.0044 mL	5.0221 mL	10.0442 mL
10 mM	0.5022 mL	2.511 mL	5.0221 mL
50 mM	0.1004 mL	0.5022 mL	1.0044 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481