


Tetradecyltrimethylammonium chloride

Chemical Properties

CAS No. :	4574-04-3
Formula:	C ₁₇ H ₃₈ ClN
Molecular Weight:	291.94
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	Tetradecyltrimethylammonium chloride is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4254 mL	17.1268 mL	34.2536 mL
5 mM	0.6851 mL	3.4254 mL	6.8507 mL
10 mM	0.3425 mL	1.7127 mL	3.4254 mL
50 mM	0.0685 mL	0.3425 mL	0.6851 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481