


[Rh(Cp-)Cl₂]₂

Chemical Properties

CAS No. :	12354-85-7
Formula:	C ₂₀ H ₃₀ Cl ₄ Rh ₂
Molecular Weight:	618.08
Storage:	Keep away from direct sunlight, Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	[Rh(Cp-)Cl ₂] ₂ is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6179 mL	8.0896 mL	16.1791 mL
5 mM	0.3236 mL	1.6179 mL	3.2358 mL
10 mM	0.1618 mL	0.809 mL	1.6179 mL
50 mM	0.0324 mL	0.1618 mL	0.3236 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481