


Tricholine citrate (65% in water)

Chemical Properties

CAS No. :	546-63-4
Formula:	C ₁₁ H ₁₉ N ₃ O ₈
Molecular Weight:	293.27
Storage:	-20°C for 2 years, 4°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tricholine citrate (65% in water) is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4098 mL	17.0491 mL	34.0983 mL
5 mM	0.682 mL	3.4098 mL	6.8197 mL
10 mM	0.341 mL	1.7049 mL	3.4098 mL
50 mM	0.0682 mL	0.341 mL	0.682 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481