

## 2-(Pyridin-2-yl)acetic acid hydrochloride

## Chemical Properties

CAS No. :	16179-97-8
Formula:	C7H8ClNO2
Molecular Weight:	173.597
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	2-(Pyridin-2-yl)acetic acid hydrochloride is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7604 mL	28.8018 mL	57.6037 mL
5 mM	1.1521 mL	5.7604 mL	11.5207 mL
10 mM	0.576 mL	2.8802 mL	5.7604 mL
50 mM	0.1152 mL	0.576 mL	1.1521 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481