


## Potassium pyrosulfate

## Chemical Properties

CAS No. :	7790-62-7
Formula:	K <sub>2</sub> O <sub>7</sub> S <sub>2</sub>
Molecular Weight:	254.32
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Potassium pyrosulfate (Pyrosulfuric acid, dipotassium salt) is a biochemical reagent that serves as a biomaterial or organic compound in life science research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9321 mL	19.6603 mL	39.3205 mL
5 mM	0.7864 mL	3.9321 mL	7.8641 mL
10 mM	0.3932 mL	1.966 mL	3.9321 mL
50 mM	0.0786 mL	0.3932 mL	0.7864 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481