


Sodium methanesulfonate

Chemical Properties

CAS No. :	2386-57-4
Formula:	CH ₃ NaO ₃ S
Molecular Weight:	118.09
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Sodium methanesulfonate [Sodium; methanesulfonate] is usable as a biomaterial or organic compound relevant to life science research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.4681 mL	42.3406 mL	84.6812 mL
5 mM	1.6936 mL	8.4681 mL	16.9362 mL
10 mM	0.8468 mL	4.2341 mL	8.4681 mL
50 mM	0.1694 mL	0.8468 mL	1.6936 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481