

1-Propanesulfonic acid sodium hydrate

Chemical Properties

CAS No. : 304672-01-3


Formula: C₃H₉NaO₄S

Molecular Weight: 164.16

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Propanesulfonic acid (sodium hydrate) is a versatile biochemical reagent utilized in various applications such as synthetic reagents and laboratory buffers. It functions effectively as a biomaterial or organic compound in life science research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0916 mL	30.4581 mL	60.9162 mL
5 mM	1.2183 mL	6.0916 mL	12.1832 mL
10 mM	0.6092 mL	3.0458 mL	6.0916 mL
50 mM	0.1218 mL	0.6092 mL	1.2183 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481