

Recombinant Humanized Type III Collagen 28.6kDa

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Recombinant Humanized Type III Collagen 28.6kDa is a synthesized form of humanized Type III collagen with a molecular weight of 28.6kDa. Known as rhCOLIII, it serves multiple biological functions, such as enhancing skin extracellular matrix regeneration and improving the cellular microenvironment. It inhibits the proliferation, migration, and invasion of breast cancer cells. Type III collagen is involved in cell adhesion, migration, proliferation, and differentiation through interactions with integrins.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481