


## TMPyP tetrachloride

## Chemical Properties

CAS No. :	92739-63-4
Formula:	C <sub>44</sub> H <sub>38</sub> Cl <sub>4</sub> N <sub>8</sub>
Molecular Weight:	820.64
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



## Biological Description

Description	TMPyP tetrachloride is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2186 mL	6.0928 mL	12.1856 mL
5 mM	0.2437 mL	1.2186 mL	2.4371 mL
10 mM	0.1219 mL	0.6093 mL	1.2186 mL
50 mM	0.0244 mL	0.1219 mL	0.2437 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481