

Tetraheptylammonium iodide


Chemical Properties

CAS No. : 3535-83-9

Formula: C₂₈H₆₀IN

Molecular Weight: 537.69

Storage: Keep away from moisture, Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tetraheptylammonium iodide is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8598 mL	9.299 mL	18.5981 mL
5 mM	0.372 mL	1.8598 mL	3.7196 mL
10 mM	0.186 mL	0.9299 mL	1.8598 mL
50 mM	0.0372 mL	0.186 mL	0.372 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481