

3'-O-Azidomethyl-dATP sodium solution (100mM)

Chemical Properties

CAS No. :

Formula: C₁₂H₁₉N₈O₁₂P₃.xNa

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3'-O-Azidomethyl-dATP sodium solution (100mM) is a labeled and modified deoxynucleoside triphosphate (dNTP) that acts as chain-elongating inhibitor of DNA polymerase for DNA sequencing.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481