


## Withaphysalin S

## Chemical Properties

CAS No. :	949172-13-8
Formula:	C <sub>29</sub> H <sub>40</sub> O <sub>8</sub>
Molecular Weight:	516.62
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Withaphysalin S is a natural product for research related to life sciences. The catalog number is TN5319 and the CAS number is 949172-13-8.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9357 mL	9.6783 mL	19.3566 mL
5 mM	0.3871 mL	1.9357 mL	3.8713 mL
10 mM	0.1936 mL	0.9678 mL	1.9357 mL
50 mM	0.0387 mL	0.1936 mL	0.3871 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481