

2,4-Dihydroxy-6-methoxy-3,5-diprenylacetophenone

Chemical Properties


CAS No. :

Formula: C19H26O4

Molecular Weight: 318.413

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,4-Dihydroxy-6-methoxy-3,5-diprenylacetophenone is a useful organic compound for research related to life sciences and the catalog number is T126162.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1406 mL	15.703 mL	31.406 mL
5 mM	0.6281 mL	3.1406 mL	6.2812 mL
10 mM	0.3141 mL	1.5703 mL	3.1406 mL
50 mM	0.0628 mL	0.3141 mL	0.6281 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481