


Zearalanone

Chemical Properties

CAS No. :	5975-78-0
Formula:	C ₁₈ H ₂₄ O ₅
Molecular Weight:	320.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Zearalenone is a nonsteroidal estrogenic mycotoxin produced by Fusarium species that colonizes various grains.
Targets(IC50)	Estrogen Receptor/ERR,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1213 mL	15.6065 mL	31.2129 mL
5 mM	0.6243 mL	3.1213 mL	6.2426 mL
10 mM	0.3121 mL	1.5606 mL	3.1213 mL
50 mM	0.0624 mL	0.3121 mL	0.6243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Zöllner P, et al. Triacsin Concentration levels of zearalenone and its metabolites in urine, muscle tissue, and liversamples of pigs fed with mycotoxin-contaminated oats. J Agric Food Chem. 2002 Apr 24;50(9):2494-501.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481