

## Mal-amido-PEG16-acid

## Chemical Properties

CAS No. :

Formula: C42H76N2O21

Molecular Weight: 945.05



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Mal-amido-PEG16-acid is a PEG-based linker used in PROTAC molecules, joining two essential ligands to facilitate selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker, one targeting an E3 ubiquitin ligase and the other the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0581 mL	5.2907 mL	10.5815 mL
5 mM	0.2116 mL	1.0581 mL	2.1163 mL
10 mM	0.1058 mL	0.5291 mL	1.0581 mL
50 mM	0.0212 mL	0.1058 mL	0.2116 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481