

Poly(ethylene glycol) diacrylate (MW 4000), MEHQ as inhibitor

## Chemical Properties

CAS No. :

Formula: C<sub>8</sub>H<sub>10</sub>O<sub>4</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Poly(ethylene glycol) diacrylate (MW 4000), MEHQ as inhibitor (PEGDA (MW 4000), MEHQ as inhibitor) is a PEG derivative used to form crosslinked polymers.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481