

Polyethylenimine (linear, average Mn 10000)

Chemical Properties

CAS No. :

Formula: (C₂H₅N)_n

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Polyethylenimine (linear, average Mn 10000, PDI≤1.3) [PEI (linear, average Mn 10000, PDI≤1.3)] is a biochemical reagent employed as a biomaterial or organic compound in life science research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481