

## Poly(ethylene glycol) methacrylate (MW 500)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Poly(ethylene glycol) methacrylate (MW 500) is a PEG functionalized with methacrylate ester, utilized in the preparation of composites such as Poly(ethylene glycol) methacrylate-chitosan, which serves as a delivery vehicle for ocular drug applications.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481