

HPMC (Type I,Viscosity:100mPa.s)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HPMC (Type I,Viscosity:100mPa.s) is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481