

## Poloxamer 334 (P104)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Poloxamer 334 P104 is a block copolymer of polyoxyethylene and polyoxypropylene with an average molecular weight of 5900. It forms a thermoreversible hydrogel and is utilized as a food additive. Additionally, Poloxamer 334 P104 serves as a drug delivery vehicle in cosmetics, pharmaceuticals, and tissue engineering.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481