

Polyacrylamide (Nonionic, Mw 5-6 million)

Chemical Properties

CAS No. :

Formula: C₃H₅NO)_x

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Polyacrylamide, Nonionic, Mw 5-6 million (Acrylamide polymer, Nonionic, Mw 5-6 million), is a versatile polymer known for its moisturizing properties. This compound can be utilized as a biomaterial or organic compound in life science research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481