

mPEG-Mal (MW 350)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

mPEG-Mal (MW 350) is a polyethylene glycol (PEG) derivative utilized for thiol PEGylation of protein molecules. The maleimide group (-Mal) degrades in aqueous environments and is employed in drug delivery research.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481