

Bovine Serum Albumin (Biotinylated)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bovine Serum Albumin (Biotinylated) is a biotin-conjugated bovine serum albumin that selectively binds to various vascular subgroups within lymphoid tissues. Serum albumin is a versatile protein with remarkable ligand-binding capabilities, making it an essential transport molecule for a range of metabolites, drugs, nutrients, metals, and other molecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481