

Dextran sulfate sodium salt (MW 9000-20000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dextran sulfate sodium salt (MW 9000-20000) (DSS (MW 9000-20000)) is a negatively charged sulfated polysaccharide. It exhibits antiviral, antibacterial, anti-inflammatory, antifibrotic, and wound healing properties. Dextran sulfate sodium salt can be used as an additive in cell culture media to prevent cell aggregation and also serves as a gel-forming agent in cosmetics.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481