

Dextran sulfate, high sulfate (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dextran sulfate, high sulfate (MW 5000), is a sulfated polysaccharide that functions as a biopolymer. It exhibits antiviral, antibacterial, anti-inflammatory, antifibrotic, and wound-healing properties. Additionally, dextran sulfate can be utilized as an additive in cell culture media to prevent cell aggregation and serves as a gelling agent in cosmetics.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481