

Dextran sulfate sodium salt (MW 6500-10000)

Chemical Properties

CAS No. : 9011-18-1

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dextran sulfate sodium salt (MW 6500-10000) (DSS (MW 6500-10000)) is a negatively charged sulfated polysaccharide. It possesses antiviral, antibacterial, anti-inflammatory, antifibrotic, and wound-healing properties. Additionally, it can be used as an additive in cell culture media to prevent cell aggregation and as a gel-forming agent in cosmetics.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481