

Methoxypolyethylene glycol maleimide (MW 750)

Chemical Properties

CAS No. :

Formula: (C₂H₄O)_nC₇H₉N₂O₉

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Methoxypolyethylene glycol maleimide (MW 750) is a reagent used in biochemical reactions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481