

DSPE-PEG-FITC (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-FITC, with a molecular weight of 3400, is a PEG lipid conjugated with fluorescein, suitable for preparing liposomes as drug carriers for targeted delivery. FITC is a green dye with a peak absorption at 494 nm and maximum emission at 520 nm, ideal for staining biological samples or nanoparticles.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481