


(2-Methoxyethoxymethyl)triethylammonium chloride

Chemical Properties

CAS No. : 60043-43-8
 Formula: C₁₀H₂₄ClNO₂
 Molecular Weight: 225.76
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2-Methoxyethoxymethyl)triethylammonium chloride is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4295 mL	22.1474 mL	44.2948 mL
5 mM	0.8859 mL	4.4295 mL	8.859 mL
10 mM	0.4429 mL	2.2147 mL	4.4295 mL
50 mM	0.0886 mL	0.4429 mL	0.8859 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481