


3-(Trifluoromethyl)phenyltrimethylammonium bromide

Chemical Properties

CAS No. :	262608-95-7
Formula:	C10H13BrF3N
Molecular Weight:	284.12
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-(Trifluoromethyl)phenyltrimethylammonium bromide is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5196 mL	17.5982 mL	35.1964 mL
5 mM	0.7039 mL	3.5196 mL	7.0393 mL
10 mM	0.352 mL	1.7598 mL	3.5196 mL
50 mM	0.0704 mL	0.352 mL	0.7039 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481