

TRIS maleate

Chemical Properties



CAS No. : 72200-76-1

Formula: C₈H₁₅NO₇

Molecular Weight: 237.21

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TRIS (maleate) is utilized in various growth media and tissue preparations, as well as serving as a diluent or wash buffer.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2157 mL	21.0784 mL	42.1567 mL
5 mM	0.8431 mL	4.2157 mL	8.4313 mL
10 mM	0.4216 mL	2.1078 mL	4.2157 mL
50 mM	0.0843 mL	0.4216 mL	0.8431 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481