

mPEG-Mal (MW 5000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store under nitrogen  
Store at 4°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

mPEG-Mal (MW 5000) is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481