

Polyvinyl alcohol (Mw 9000-10000)

## Chemical Properties

CAS No. :

Formula: (C<sub>2</sub>H<sub>4</sub>O)<sub>x</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Store at RT

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

Polyvinyl alcohol (Mw 9000-10000, 80% hydrolyzed), also known as PVA, is a water-soluble vinyl polymer that functions as a nonionic surfactant. Additionally, it serves as a biodegradable polymer suitable for use in adhesives, coatings, textiles, ceramics, and cosmetics.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481