

Poloxamer 188 (F68)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Poloxamer 188 is a nonionic linear copolymer known for its surfactant properties. In various tissue damage models, Poloxamer 188 demonstrates antithrombotic, anti-inflammatory, and cytoprotective activities.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481