

PEG10000

## Chemical Properties

CAS No. :

Formula:

Molecular Weight: 10000

Storage:

Store at RT

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PEG10000 (Polyethylene glycol 10000) is utilized as a solubilizing agent and can be employed as a biomaterial or organic compound in life sciences research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.100 mL	0.500 mL	1.000 mL
5 mM	0.020 mL	0.100 mL	0.200 mL
10 mM	0.010 mL	0.050 mL	0.100 mL
50 mM	0.002 mL	0.010 mL	0.020 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481