

PEG20000

## Chemical Properties

CAS No. :

Formula:

Molecular Weight: 20000

Keep away from direct sunlight

Storage:

Store at RT

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PEG20000 (Polyethylene glycol 20000) can be used as a solubilizer. It serves as a biomaterial or organic compound relevant to life science research and is biodegradable by bacteria.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.050 mL	0.250 mL	0.500 mL
5 mM	0.010 mL	0.050 mL	0.100 mL
10 mM	0.005 mL	0.025 mL	0.050 mL
50 mM	0.001 mL	0.005 mL	0.010 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481