

Polyvinylpyrrolidone (average Mw ~ 360000)

## Chemical Properties

CAS No. :

Formula: (C<sub>6</sub>H<sub>9</sub>NO)<sub>n</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Polyvinylpyrrolidone (average Mw ~ 360000) is a versatile synthetic polymer with an average molecular weight of 360000. It is extensively utilized in the synthesis of nanoparticles.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481